


Dosimetry Check™


3D | EPID Based | Pre/In-Vivo/Exit Dosimetry | Automated QA

EPID Based IMRT and Exit Dosimetry

Collapsed Cone Convolution Superposition Algorithm



- EPID Based Measurement
- Monitor Unit Verification
- Independent Algorithm: Uses FDA 510(k) cleared collapsed cone calculation engine
- No Diode Array Setup Required
- No Phantom Setup Required



Dosimetry Check Uses Fluence to Calculate Dose

Sales Support Inside USA: 866-592-1343 option 1 • Outside USA: +1-903-207-4298 option 1 • Technical Support: option 2
Email: sales@lifelinesoftware.com Website: www.lifelinesoftware.com

*Dosimetry Check users with active contract will be automatically upgraded to the *RadCalc 3D platform.

*RadCalc 3D Platform Is Not Available For Sale Within The United States – Work in Progress

MKT-013 REVISION D

