

Dosimetry Check™

3D | EPID Based | Pre/In-Vivo/Exit Dosimetry | Automated QA


Providing Essential QA Tools TomoTherapy DQA

Collapsed Cone Convolution Superposition Algorithm



- Detector Based Measurement
- Monitor Unit Verification
- Independent Algorithm
- No Diode Array Setup Required
- No Phantom Setup Required

Models the TomoEdge Dynamic Jaw Movement Providing Essential QA Tools



Sales Support Inside USA: 866-592-1343 option 1 • Outside USA: +1-903-207-4298 option 1 • Technical Support: option 2
Email: sales@lifelinesoftware.com Website: www.lifelinesoftware.com

*Dosimetry Check users with active contract will be automatically upgraded to the *RadCalc 3D platform.

*RadCalc 3D Platform Is Not Available For Sale Within The United States – Work in Progress

MKT-002 REVISION D